

Proper care for TCS appliances: To keep your appliance looking and feeling like new, please follow these simple recommendations:

1. Continue to brush your natural teeth and gums as directed by your dentist.
2. Clean your appliance daily. Rinse under warm running water for approximately 10-30 seconds to remove any food particles.
3. Brush your appliance with a soft denture brush using water only. Rinse again.
4. Place your appliance in denture cleaner for at least 15 minutes every day. We recommend tcs® Concentrated Dental Appliance Cleaner.
5. When possible, keep your appliance in water or a denture solution when you are not wearing it. DO NOT use bleach, toothpaste, or similar cleaning products, these may scratch or discolor the surface of your appliance.

- A)**  **TCS Dental Appliance Cleaner (Item #4505-01)**  
Deep cleans, removes stains, and kills bacteria and odors. Used daily, it will keep your appliance looking and feeling like new!  
Box of 10; approx. a 10 month supply.
- B)**  **TCS Sonic Cleaner (Item #4013-01)**  
Obtain professional cleaning results with the aid of sonic vibrations. Must be used in combination with a dental appliance cleaner. Also, ideal to use as a storage container. (2AA batteries included).
- C)**  **TCS Cool Mint Spray (Item #4504/PY-01)**  
Mint flavored spray used to clean and refresh your dental appliance on-the-go. Yields approx. 200 misty sprays.
- D)**  **Dental Appliance Brush (Item #4506FR-01)**  
Soft bristle brush to thoroughly clean dental appliances without scratching the surface.

**Order options:**  
**Call:** 866-426-2970  
**Online:** www.tcsfresh.com  
**Mail:** Send form and payment  
**to:** TCS, Inc. 2619 Lime Ave. Signal Hill, CA 90755  
**Ship to USA only**

- a) TCS Dental Appliance Cleaner **\$ 30.00**
- b) Sonic Cleaner **\$ 18.95**
- c) Cool Mint Spray **\$ 12.00**
- d) Dental Appliance Brush **\$ 3.00**

Save with packages

- A + B **\$ 43.95** **save \$ 5.00**
- A + B + C **\$ 52.95** **save \$ 8.00**
- A + B + C + D **\$ 54.95** **save \$ 9.00**

TOTAL: \_\_\_\_\_ FREE SHIPPING

- check enclosed     amex     master card / visa     discover

card#	exp. date	CVV
name	phone #	
street	apt #	
city	state	zipcode
email		

## FAQ's

### How long will the procedure take?

Sometimes a partial or denture can be delivered in as little as two visits, but usually the average is 4 visits. Approximately two to six weeks.

### Will the process hurt?

TCS restorations are non-invasive, so there is minimal discomfort. For partial dentures, sometimes preparations to the abutment teeth are necessary, but there is little discomfort associated with this procedure.

### How long will my denture last?

TCS restorations are designed to last a lifetime. Should modifications be needed as time goes by due to natural changes in your mouth, your partial or denture can usually be modified rather than replaced.

### Can I wear my denture at night?

Some designs of dentures may not be recommended for overnight use. Ask your dentist.

### Will I have difficulty chewing?

Some foods may never be easy to chew with dentures, but a properly designed and fitted denture combined with TCS material's durability and superior adaptation will give you the best chewing experience of any denture material available.

### Will the denture change how I speak?

Sometimes it is difficult to speak clearly with missing teeth, therefore, the denture may improve your speech clarity. If you find it difficult to pronounce certain words, do not worry, over time you will become accustomed to speaking properly with your new appliance.

### Can I make minor adjustments to my denture?

We strongly recommend that you do not alter you appliance in any way. If adjustments are needed, please contact your dentist.

### How long will it take me to adapt to my new denture(s)?

Partial dentures are normally easier to become accustomed to because they attach to your natural teeth. Full dentures require muscles to work in new ways in order to help stabilize and manipulate them. Denture wearers can feel completely comfortable with just a little practice.

### Should I consider implants?

Dental implants are a possible alternative to dentures. Implants are surgically integrated into your jaw bone to hold a replacement tooth or bridge in place. Implants and bridges may resemble the "feel" of real teeth, however they are significantly more expensive than dentures and the process is invasive. Not all patients are good candidates for implants, so be sure to talk to your dentist about which treatment option is best for you.



# Removable Dentures

Natural, Comfortable & Affordable

- Learn how removable dentures may be an ideal solution for missing teeth
- Tips for Denture Care and Maintenance
- FAQ



www.tcsdentalinc.com

8010-A1 Rev E

# Ask your dentist if TCS products are right for you.



Full Dentures



1-2 Tooth Appliances



Partial Dentures

Adults are living longer and healthier lives, it is not unusual that people will lose one or more teeth in their lifetime, even more so as life expectancy continues to increase. Luckily there are plenty of options for tooth replacement, you do not have to worry about a denture compromising your smile. For many patients that are concerned about the unreliability of temporary appliances such as flippers, revealing partials with metal clasps, or expensive options like implants, that require invasive surgery and long healing periods, TCS removable dentures may be the ideal choice. When it comes to deciding on a replacement, many things are considered – i.e. cost, procedure, aesthetics, etc. These dental care decisions should always be made in consultation with your dentist.

Partial dentures fill in the gaps created by missing teeth and help prevent your natural teeth from shifting. You will be happy to hear that a well-fitted partial denture can provide wonderful aesthetics and functionality. TCS flexible partials are usually fabricated without unsightly metal clasps that will compromise your smile. In some cases, metal may be incorporated, this is a functionality decision made by your dental professional.

TCS flexible partials offer an ideal degree of flexibility which allows the partial(s) to flex with the natural contours of the mouth for maximum comfort. TCS offers two base materials for partial dentures: **iFlex** and **Unbreakable**. Ask your dentist which is the right choice for you.

Full dentures replace all missing teeth of either the upper arch, lower arch or both. Full dentures function between the muscle forces of the lips, cheeks and tongue. Wearing a well-fitted denture is not only essential to your health, it will also rejuvenate your appearance. **Karadent** is a technologically advanced resin used for the fabrication of full dentures. Ill fitting dentures can slip and move around causing discomfort and embarrassment. **Karadent** is processed using high pressure injection for a precise fit. **Karadent** is monomer-free, making it an ideal choice for patients sensitive to traditional acrylic dentures.

Unlike common belief, dentures are not just for elderly patients. Patients of any age may lose some or all of their teeth and need an affordable, aesthetic solution. You may be surprised at the amount of people under the age of 35 that wear dentures as a permanent solution.



A 1.5 minute animation video that demonstrates how a typical TCS partial denture compares to a conventional metal and acrylic partial.



TCS is featured in "Improving Lives through Effective Dental Solutions," a 5 minute Health Brief hosted by NFL star, Terry Bradshaw.

Scan to view videos about TCS products or visit:  
[www.tcsdentalinc.com](http://www.tcsdentalinc.com)



## Tiffany's Story

Tiffany lost most of her teeth in an automobile accident at the age of 12. After many years of wearing an uncomfortable metal and acrylic partial denture that would constantly break, Tiffany was finally prescribed an Unbreakable by tcs®. Tiffany is excited about finally wearing something she can count on, looks great, and feels even better.



TCS Dental Appliances are comfortable, easy to care for, and best of all no one can tell you are wearing one.

See reverse for tips on denture care

Tooth loss can have an emotional impact on some people, please talk to your dentist about emotions you are feeling about tooth loss.

Excellent dental oral hygiene, proper maintenance of your restoration and regular dental checkups are essential for your health.