

Ideal Choice for Removable Dental Appliances



TIP:

Hot Water Treatment for Flexible Partial

The Hot Water Treatment is very important as it allows for a smooth initial insertion and a good adaptation with your patient's natural tissue. (This is only necessary for 1st insertion).

1. Examine the partial carefully on and off the master model and check the path of insertion.
2. Submerge the TCS partial denture in very hot tap water for approximately one minute.
3. When the partial cools to a tolerable temperature, gently insert it into the patients mouth. You may choose to spray the partial with TCS Cool Mint Spray prior to initial insertion.



“ My partial is absolutely undetectable, it is thin and comfortable, and it fits like a glove. I feel completely confident with my TCS partial. ”

-Mike Hall | 52 years old | NC



Sample Models

Unbreakable™ partial 4011-01 ea.
iFlex™ partial 4009-01 ea.
Karadent™ full denture 4016-01 ea.



Green Mounted Stones

Ideal for chairside adjustments.
4501-01 box of 100



Brown Rubber Points

Ideal for chairside adjustments. Use after green mounted stone, when adjusting Unbreakable™ and/or iFlex™.
Small 4511-01S ea.
Large 4510-01L ea.



Carbide Bur

Specifically designed to cut through thermoplastic materials. These burs are long lasting with high cutting efficiency.
50C79-01 ea.



Carving Knife & Replacement Blades

Ideal for removing fibers that sometimes appear on thermoplastic restorations during adjustments. Replacement blades are specific to this knife.
Knife 4201-01 ea. | Blades 4202-01 pack of 12



Mini Rag Buff

Recommended for touch-up polishing of TCS appliances. Use with blue hi-shine.
4221-01 ea.



Mini Fuzz Buff

Recommended for touch-up polishing of TCS appliances.
4222-01 ea.



Blue Hi-Shine

Recommended for touch-up polishing of TCS appliances. Use with the mini rag buff.
4215-01 (6.7 oz) ea.

Dr. Adjustment Kit

Includes: 1 green mounted stone, 1 lg brown rubber point, 1 sm brown rubber point, 1 carbide bur, 1 carving knife, sample piece of blue hi-shine paste, 5 mini rag buffs, 5 mini fuzz buffs.
4014-02

ORDER OPTIONS:

866-426-2970 | www.tcsdentalinc.com

TCS Dental Appliance Cleaner

Deep cleans, removes stains, and kills bacteria and odors. Used daily, it will keep appliances looking and feeling like new.

4505-01 **Box of 10** (4 week supply packets)
* Ideal for patient purchase. 10-month supply.
4505FR-01 **Box of 70** (4 week supply packets)
* Ideal to give as samples
4505FR-SMP **Box of 75** (1 week supply packets)
* Ideal to give as samples



TCS Sonic Cleaner

Patients can obtain professional cleaning results at home with the aid of sonic vibrations. Must be used in combination with dental appliance cleaner. Also, ideal to use as a storage container. (2AA batteries included)
4013-01



TCS Cool Mint Spray

Mint flavored spray used to refresh dental appliances on-the-go. Yields approx. 200 misty sprays.
4504/PY-01



TCS Dental Appliance Brush

Soft bristle brush to thoroughly clean dental appliances without scratching the surface.
4506FR-01