

# Exceptionally strong and reliable. Over 10 years IPS e.max.

IPS e.max Press is based on more than ten years of experience, profound expertise and innovative strength. We are driven by our passion to produce esthetic teeth and the pleasure of seeing happy and smiling patients. Various long-term studies confirm the safety and impressive reliability of the material. As a result, patients can rely on restorations that will last for many years.



Starting situation



After seating the restorations



After 10 years *in situ*

Dr Sidney Kina / José C. Romanini, Brazil

## Maximum dependability

IPS e.max Press is characterized by its high strength. Quality tests of more than ten years show that IPS e.max Press demonstrates a mean biaxial strength of 470 MPa<sup>4</sup>. This value by far surpasses the required standard of 100 MPa for crowns placed with the adhesive technique and 300 MPa for conventionally placed crowns.

”

I've been using IPS e.max Press for more than ten years and I can assure you that the material is one of the best products that I've had the privilege of using.

”



August Bruguera, Spain

**97.5%**  
survival rate<sup>1</sup>

**10** years of  
clinical evidence

More than **100**  
million restorations<sup>2</sup>

**98%**  
customer satisfaction<sup>3</sup>

Strength of  
**470 MPa**<sup>4</sup>

<sup>1</sup>IPS e.max Scientific Report Vol. 02/2001-2013

<sup>2</sup>IPS e.max System, based on sales figures


<sup>3</sup>Corporate Marketing Insight Ivoclar Vivadent, Schaan, Liechtenstein

<sup>4</sup>Mean biaxial strength,

R&D Ivoclar Vivadent, Schaan, Liechtenstein

## Fixed prosthetics

IPS e.max® forms a part of the "Fixed Prosthetics" product category. The products of this category cover the procedure involved in the fabrication of fixed prosthetic restorations – from temporization to restoration care. The products are optimally coordinated with each other and enable successful processing and application.



### THESE ARE FURTHER PRODUCTS OF THIS CATEGORY:

#### Programat®

Press and ceramic furnaces for demanding requirements



Packed with proven technology and advanced innovations

- Outstanding press and firing results
- Ideally coordinated with the ceramic materials of Ivoclar Vivadent
- Easy operation

#### Variolink® Esthetic

The esthetic luting composite



The luting composite for exceptional esthetics and user-friendly processing

- Balanced and concise Effect shade system
- Excellent shade stability due to amine-free composition
- Easy, controlled excess removal

Would you like to know more about the products of the "Fixed Prosthetics" category?

Simply get in touch with your contact person at Ivoclar Vivadent or visit [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com) for more information.

Production and distribution  
**Ivoclar Vivadent AG**  
Bendererstr. 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Descriptions and data constitute no warranty of attributes.  
Printed in Germany  
© Ivoclar Vivadent AG, Schaan/Liechtenstein  
690520/en/2017-02-27

**ivoclar**  
**vivadent**  
passion vision innovation